

NORHOF Liquid Nitrogen Dewars

Aluminum Dewars for Liquid Nitrogen storage

- 💧 Aluminum Dewars designed for storage and transport of Liquid Nitrogen
- 💧 All units are CE certified
- 💧 Static Dewar type (no pressure build up)
- 💧 Small Neck with NW50 flange



Technical Specifications

| Specifications | 900-000-012 | 900-000-025 | 900-000-035 | 900-000-050 | 900-000-100 |
|-----------------------------|-------------|-------------|-------------|-------------|-------------|
| Capacity (l) | 12 | 25 | 35 | 50 | 100 |
| Neck Diameter (mm) | 51 | 51 | 51 | 51 | 51 |
| Total Height (mm) | 571 | 681 | 602 | 717 | 1035 |
| Inner Height (mm) | 533 | 612 | 533 | 640 | 922 |
| Outer diameter (mm) | 300 | 400 | 460 | 460 | 500 |
| Empty weight (kg) | 7 | 12 | 13 | 15 | 32 |
| Full weight (kg) | 16 | 33 | 40 | 57 | 113 |
| Static holding time (days) | 84 | 104 | 112 | 122 | 164 |
| NW50 flange | Yes | Yes | Yes | Yes | Yes |
| Max. working pressure (bar) | 0,5 | 0,5 | 0,5 | 0,5 | 0,5 |
| Dewar trolley option | 900-400-012 | 900-400-025 | 900-400-000 | 900-400-000 | 900-400-000 |

Dewar Trolley Option:

